

# VignettePTRViewer

## Premiliminary work

This vignette shows how to use chemosensR package to read .txt files, calculate AUC and return statistical analyses. This command updates the chemosensR package

```
library(PTRMSR)
library(devtools)
library(MSnbase)
library(reshape2)
library(pheatmap)
library(ggplot2)

library(ellipse)
install_github("https://github.com/ChemoSens/ChemoSensPrivate/CSUtils")
install_github("https://github.com/ChemoSens/PTRMSR")
```

Then, chemosensR should be loaded.

```
library(PTRMSR)
```

## Meta-data

All the data files (.txt from ptrViewer) and the metadata file (.csv) should be in a single repository (wd)

```
wd="C:/INRA/Data/Donnees Cantin/Cantin-DTS-PTRviewer"
setwd(wd)
listFiles=list.files(pattern="*.txt")[-1]
metaData2=read.table("metaData2.csv",sep=";",header=T)
head(metaData2[, -c(2:3)])
```

## Analysis of PTR-Viewer files

### Selection of relevant ions only

This function returns the names of ions whose maximal intensity is 3 times higher than the maximal intensity during the noise period. It returns a list containing (i) a list containing for each file the ratio maximal intensity during the tasting/maximal intensity during the noise period, (ii) a vector (intersection) containing the ions which are significant in all files (iii) a vector (union) containing the ions which are significant in at least one file

```
setwd(wd)
referenceBreath="m69.06906..69.06906...Conc."
sigIons=ptrvListSignificantSNRIons(listFiles=listFiles[1:4],
                                   metaData=metaData2,noisePeriod=c(0,25))
ionSigUnique=sigIons$union
```

To obtain further help for this (or any) function, enter:

```
?ptrvListSignificantSNRIons
```

## AUC calculations

This function allows statistics to be calculated for each file after breathing correction

The results can be saved into csv files

```
res_auc$listRes
write.table(file="auc.csv",sep=";",res_auc$listRes,row.names=F)
```

The stat option can be used to select the statistic.

```
res_tmax=ptrvListIntensityByTime(listFiles=listFiles,metaData=metaData2,ions=ionSigUnique,stat="tmax")
```

## Statistical analysis

In this part, the evaluations can be normalized (a total abundance of 1 for each evaluation) or logged. ###  
PCA of observations

```
respca=ptrvListPCA(ptrvList$totalIntensity,dataType="productMeans",log=FALSE)
# rajouter projections individuelles plutot qu'ellipses
plotPCAagg(respca,type="ind")
plotPCAagg(respca,type="var",text=FALSE)
p=plotPCAagg(respca,type="var",text=FALSE)
p
library(plotly)
ggplotly(p)

p=plotPCAagg(respca,type="cor1",n=80)
plotPCAagg(respca,type="cor2")

respca=ptrvListPCA(ptrvList$totalIntensity,dataType="raw")
p=plotPCAagg(respca,type="ind",text=FALSE)
ggplotly(p)
plotPCAagg(respca,type="var",text=FALSE)
```

## Anova results

```
resanova=ptrvListAnova(ptrvList$totalIntensity[ptrvList$totalIntensity[, "ion"]=="m101.09"],normalizeBy)
```

## Heatmaps

```
ptrvListHeatmap(df=ptrvList$totalIntensity,formula=as.formula(product+subject+rep~ion),
  fun.aggregate="mean",clusterRows=T,clusterCols=T)
ptrvListHeatmap(df=ptrvList$totalIntensity,formula=as.formula(product+subject+rep~ion),
  fun.aggregate="mean",clusterRows=T,clusterCols=T)
ptrvListHeatmap(df=ptrvList$totalIntensity,formula=as.formula(product+subject+rep~ion),
  fun.aggregate="max",clusterRows=T,clusterCols=T)
ptrvListHeatmap(df=ptrvList$totalIntensity,formula=as.formula(product+subject~ion),
  fun.aggregate="max",clusterRows=T,clusterCols=T,showRownames=T)
ptrvListHeatmap(df=ptrvList$totalIntensity,formula=as.formula(product+subject~ion),
  fun.aggregate="mean",normalization="row",clusterRows=T,clusterCols=T,showRownames=T)
ptrvListHeatmap(df=ptrvList$totalIntensity,formula=as.formula(product+subject~ion),
  fun.aggregate="mean",normalization="row",clusterRows=T,clusterCols=T,showRownames=T,ann
```

## Outputs for one given file

It could be interesting to see the raw data of one file. In this purpose, the function `ptrvReport` returns several plots which could help to interpretation.

```
setwd(wd)
file=listFiles[1]
dataset=read.table(file=file,header=TRUE,sep="\t")
report=ptrvReport(dataset,selectIons="evolving",
                  listIons=ionSigUnique[1:3],
                  referenceBreath=referenceBreath,
                  methodDetectStart="startPeakProportion",
                  noisePeriodIBT=c(0,30),noisePeriodSig=c(0,30),
                  noisePeriodDS=c(0,30),
                  proportionOfMax=0.3,halfWindowSize=12,maxPeaks=30)

names(report$gg)
report$gg$p_breath
```

Regarding the breathing cycle detection, these two outputs gives the results of cycle limits and the smooth data for breathing (useful to adjust the parameters `halfWindowSize` and `maxPeaks`)

```
plot(report$gg$p_breath$p_cyclelimits)
plot(report$gg$p_breath$p_smoothbreath)
```

Regarding the ions distributions

```
plot(report$gg$p_curves$p_raw)
plot(report$gg$p_curves$p_cycle)
```

Regarding the breathing ion distribution

```
plot(report$gg$p_curves$p_breath_raw)
plot(report$gg$p_curves$p_breath_cycle)
```

Do not hesitate to use `plotly` for interactive graphs

```
library(plotly)
ggplotly(report$gg$p_curves$p_breath_raw)
```

## Appendix

### Help for filling meta-data

```
ptrvCreateMetaFile(wd,subject=c(9,12),product=NULL,replicate=NULL,sep=";")
metaData=read.csv("metaData.csv",sep=";",header=TRUE,dec=",")
metaData[, "resp"]=referenceBreath
```

###[OPTIONAL] Completing meta-data file with the tasting time

```
# courbes + lissage
t0=tn=suj=ordre_produit=product=repet=finTS=rep(NA,nrow(metaData))
names(t0)=listFiles
for(i in 1:nrow(metaData))
{
  print("dataset")
  print(i)
  dataset=read.table(listFiles[i],sep="\t",dec=dec_vec[i],header=T)
  dataset[, "RelTime"]=as.numeric(dataset[, "RelTime"])
```

```

# dataset[,ions]=apply(dataset[,ions],2,as.numeric)
res_intensity=ptrvIntensityByTime(dataset,ions=ionSigAll
                                ,referenceBreath=referenceBreath,
                                correction = "cycle",
                                removeNoise=TRUE,breathRatio =FALSE,timeBlank=c(0,25),
                                halfWindowSize=12,maxPeaks=30,method="SuperSmoother",total=FALSE)
t0[i]=ptrvDetectStart(res=res_intensity,starts=ionSigAll, method="startPeakProportion",proportionOfI
                    noisePeriod=c(0,25),startPeriod=c(25,100))$tx
tn[i]=max(res_intensity$res[, "time"])
suj[i]=substr(listFiles[i],9,12)
ordre_produit[i]=substr(listFiles[i],14,14)
tds_raw_i=tds_raw[tds_raw[, "SubjectCode"]==suj[i]&tds_raw[, "Rang.produit"]==ordre_produit[i],]
product[i]=unique(tds_raw_i[, "ProductCode"])
repet[i]=unique(tds_raw_i[, "Replicate"])
finTS[i]=max(tds_raw_i[, "Time"])
}
metaData[, "miseEnBouche_s"]=t0;metaData[, "finPTR_s"]=tn;
metaData[, "finTimeSens_s"]=finTS
metaData[, "product"]=product;metaData[, "subject"]=suj;metaData[, "rep"]=repet
write.table(metaData,file="metaData2.csv",sep=";",row.names=F,dec=".")

```

##Installing chemosensR (to be done only one time) As a first step, chemosensR package should be installed. As chemosensR is on a private GitHub repository, this step required the use of a token which is send to the user.

```

install.packages("devtools")
library(devtools)
install_github("MahieuB/MultiResponseR")
if (!requireNamespace("BiocManager", quietly = TRUE))
  install.packages("BiocManager")
BiocManager::install("MSnbase")
BiocManager::install("rhdf5")
install.packages("multcompView")
#BiocManager::install("mzR")
install_github("https://github.com/ChemoSens/ChemoSensPrivate/ChemosensR", auth_token="d92a

# FOR SHINY
setwd(paste0(.libPaths()[1], "/chemosensR/extdata"))
then open the .R file and run

```

Then, the package has to be loaded with library function. Contrarily to the previous installation, this step of loading is required at any use of chemosensR package.