

Using PTRMSR for automatizing pre-processing

Contents

General principle	1
Loading the library	1
Loading metadata and list of files to be processed	1
AUC calculations of several PTR-Viewer files with noise and breathing correction	1
Documentations on the functions are available in the R documentation	2

This vignette shows how to use PTRMSR package to read .txt files, calculate AUC and return statistical analyses.

General principle

In this package, the functions beginning by ptrv use a single dataset as input. The functions beginning by ptrvList use names of several files to be used and a metaData file containing information about the different files. The objective of this paper is to use files produced by ptrViewer (or any other preprocessing of PTR-MS device returning files with ion, intensity and time as columns) then return curves and relevant statistics on the evaluations. Several biases can be corrected: (i) the breathing of the subject by calculating statistics by breathing cycle (use the parameter correction="cycle"), (ii) blank period during the evaluation that can be removed from the intensities (use blankCorrection = TRUE and noisePeriod=...). Some statistics are returned.

Loading the library

```
library(PTRMSR)
```

Loading metadata and list of files to be processed

All the data files (.txt from ptrViewer) and the metadata file (.csv) should be in a single repository

```
repo="./../inst/extdata"
listFiles=list.files(repo,pattern="*.txt")
metaData=read.table(paste0(repo,"/metaData.csv"),sep=";",header=T)
head(metaData)
```

##	file	rep	breathing	into	start	nothing	stop	product	subject
## 1	file1.txt	1	m69.0690	36.93347	0	NA	129.8109	B	S002
## 2	file2.txt	1	m69.0690	39.02045	0	NA	140.6618	A	S002

AUC calculations of several PTR-Viewer files with noise and breathing correction

This function allows statistics to be calculated for each file after breathing correction

The results can be saved into csv files

```
write.table(file="auc.csv",sep=";",res_auc$listRes,row.names=F)
```

The stat option can be used to select the statistic.

```
res_tmax=ptrvListIntensityByTime(listFiles=listFiles,metaData=metaData,ions=c("m31.0183", "m31.0259",  
                                                                              "m33.0324"),stat="tmax")
```

The breathing cycles can be observed with this command.

```
grid.arrange(grobs=res_tmax$cycleLimits[1:4])
```

Documentations on the functions are available in the R documentation

```
?detectCycle  
?ptrvIntensityByTime # for only one dataset,  
?ptrvIntensity # for summarizing the time of one dataset  
?ptrvListIntensityByTime # for only one dataset,  
?ptrvListIntensity # for summarizing several datasets
```